

ROOTT implant library for Exocad software

Content

ROOTT implant library - installation instruction for Exocad software	3
Digital ROOTT details and instructions of use	8
Scan posts	10
Cement gap	11
Screw channel	12
Digital analog	14
ROOTT R Scan post	15
Abutments	16
Instruction for Exocad using ROOT R	18
ROOTT M ROOTT P Scan post	X
Abutments	X
Instruction for Exocad using ROOT M ROOTT P	X
ROOTT S Scan post	X
Abutments	X
Instruction for Exocad using ROOT S	X
ROOTT C ROOT CS ROOTT B ROOTT BS Scan post and Transfer	X
Abutments	X
Instruction for Exocad using ROOTT C ROOT CS ROOTT B ROOTT BS	X
Model Creator	X



ROOTT implant library –

installation instructions for


Exocad software

Instructions will help you to locate suitable ROOTT library folders with implants for
Exocad libraries


Instructions are available in video format at
www.youtube.com




Step 1. Find and open **EXOCAD** folder. Inside the folder locate and open **exocad- ChairsideCAD-64Bit**




Step 2. Open **ChairsideCAD** folder




Step 3. Open **library** folder




Step 4. Open **Downloads** folder (where ROOTT Library folder was downloaded) in new window




Step 5. Open **ROOTT Library** and open **Roott DentalCAD** folder




Step 6. Copy all files from **Roott DentalCAD** files




Step 7. Paste in **implant** folder (EXOCAD > exocad-ChairsideCAD-64Bit > ChairsideCAD > library > implant)




Step 8. Go back to **library** folder, choose **modelcreator** folder and then **implants** folder (implant < library > modelcreator > implants)



Step 9. In the second window go back to **ROOTT Library**, open **Roott ModelCreator** folder and copy 4 located folders (Roott DentalCAD > ROOTT Library > Roott ModelCreator)



Step 10. Paste the copied folders in **implant** folder opened in the other window



 ROOTT library is now available

For question please contact: dalia.petkeviciute@trate.com

Digital ROOTT details and instruction of use



We recommend to follow these rules :

- Do not mix different product between implant categories. The only exception is ROOTT R library in which ROOTT S and ROOTT M detail assemblies are integrated.
- ROOTT R, ROOTT M, ROOTT S, ROOTT P implant systems intended for two-piece implants, ROOTT C, ROOTT CS, ROOTT B, ROOTT BS – for one piece.
- Two piece implants works in principle following implant + together with abutment , one piece implants are based by telescopic abutment.
- Library name are created from terms: Manufacturer+ Implant + abutment, material.

DISCLAIMER OF WARRANTIES.

ALL SERVICES AND TRATE MATERIALS ARE PROVIDED “AS IS” AND “AS AVAILABLE” WITHOUT WARRANTY OF ANY KIND AND TRATE HEREBY DISCLAIMS ALL WARRANTIES, WHETHER EXPRESS, IMPLIED, STATUTORY, OR OTHER, AND TRATE SPECIFICALLY DISCLAIMS ALL IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, TITLE, AND NON-INFRINGEMENT, AND ALL WARRANTIES ARISING FROM COURSE OF DEALING, USAGE, OR TRADE PRACTICE.

WITHOUT LIMITING THE FOREGOING, TRATE MAKES NO WARRANTY OF ANY KIND THAT THE SERVICES, TRATE MATERIALS, OR THE INFORMATION CONTAINED IN THE SERVICES WILL MEET USER’S OR ANY OTHER PERSON’S REQUIREMENTS, OPERATE WITHOUT INTERRUPTION, ACHIEVE ANY INTENDED RESULT, BE COMPATIBLE OR WORK WITH ANY SOFTWARE, SYSTEM, OR OTHER SERVICES, OR BE SECURE, ACCURATE, TIMELY, RELIABLE, RELEVANT, COMPLETE, FREE OF HARMFUL CODE, OR ERROR FREE.

Scan Posts

Scan posts indicate the exact position of the implant in the jaw. During the scanning process, the information about the position is transferred into digital format.


ROOTT library has intraoral and extraoral scan posts. The difference between them is size and shape which provide better performance for a particular workflow. The scheme below shows suitable scan posts for a particular implant type. Working on with ROOTT R implants, scan posts (SPCOM, SPCOMIO, SPCOMS, SPCOMIOS) can be used by assembling them with abutments M1 or MS1. Another advantage of digital ROOTT library is that transfer and telescopic abutments of ROOTT C, CS, B, BS implants also can be used as scan posts. These mentioned possibilities create wider applicability of the products.




ROOTT M ROOTT P



ROOTT S



ROOTT C ROOTT CS ROOTT B ROOTT BS




Cement gap

An empty gap between a titanium base and an inner crown surface is called a cement gap. A dimension in microns (μm) indicates an offset from the titanium base.

Cement gap size results depends on the equipment. To get the best results, it is recommended:

- to find the most suitable cement gap option for your equipment;
- to check regularly the mill in the machine;
- to comply with the requirements of the equipment manufacturers;
- to check if the setting of prosthesis and manufacturing equipment match (see picture below).



/izard Model Alignment

MODEL TYPE:

Baumann Modellsystem 2000® EVO for medium jaw

Baumann Modellsystem 2000® EVO for small jaw

Formlabs Form2 Hollow 100 μm

Formlabs Form2 Hollow 50 μm - DentalModel resin


Formlabs Form2 Solid 100 μm

Formlabs Form2 Solid 50 μm - DentalModel resin

ROOTT digital library is created for 6 cement gap varieties:

- Zirconia - 0.025, 0.035, 0.050 μm ;
- Plastic PEEK or PMMA - 0.075, 0.090, 0.110 μm .

If there is a need to create from PEEK PMMA with a wider cement gap, choose Zro section.
Metal is adapted only for ROOTT M and S Multi Unit.




Screw channel


There are two options to create an aesthetical prosthesis:

- Standard screw channel that matches the diameters of the screw and the screwdriver.
- Narrow screw channel that matches only the screwdriver. Choose this option if there is less occlusal surface or a more aesthetical result is desired.


Screw channel sizes are available for ROOTT R, ROOTT M, ROOTT P, and ROOTT S.



Angled screw channel position can be determined by scan post. The front plane of the scan post corresponds to the angle direction of the screw channel. Therefore, we recommend turning the scan post to lingual surface or palatal surface direction before scanning in order to achieve the best results.



The front plane of scan posts



Correct scan post direction highlighted green, incorrect - red

Digital analogs

Digital analogs indicate the exact position of the implant in the jaw. Therefore, if the scan post is not screwed properly, it could lead to inaccurate position of digital analog. Digital analogs could be used only with printed models.

ROOTT **R**



AND




ROOTT **P**



ANMD

ROOTT **M**



ROOTT **S**



ANMSD



ROOTT **C**



ANED

ROOTT **CS**

ROOTT **B**

ROOTT **BS**



ROOTT R

Scan posts



Scan posts indicate the exact position of the implant in the jaw. During the scanning process, the information about the position is transferred into digital format. According to scan post position, height and direction, scan post is converted to abutment and sets the precise location of analog. Get yourself acquainted with all types of ROOTT R scan posts, which are developed for the effective and precise workflow of implantologists (intraoral) and dental technicians (extraoral).

SPCO Scan post characteristic

- Extraoral scanning
- Long body allows comfortable usage and precise results of laboratory workflow
- Easily scannable
- Reusable




SCP0IO Scan post characteristic

- Intraoral scanning
- Short body allows comfortable usage and precise results of implantologists' workflow
- Easily scannable
- Reusable



SPCOM+M1/ SCPOMS+MS Scan post characteristic

- Extraoral scanning
- Long body allows comfortable usage and precise results of laboratory workflow
- Wider applicability of SPCOM/SPCOMS (Scan post for ROOTT M, ROOTT P, ROOTT S implants)
- Easily scannable
- Reusable



SPCOMIO+M1/ SPCOMIOS+MS1 Scan post characteristic

- Intraoral scanning
- Short body allows comfortable usage and precise results of implantologists' workflow
- Wider applicability of SPCOMIO, SPCOMIOS (Scan post for ROOTT M, ROOTT P, ROOTT S implants)
- Easily scannable
- Reusable



Abutments

ROOTT R digital abutments are suitable for ordinary and complicated clinical situations. There is a wide range of options for multi-unit and single crown cases.

CRE is a multi-functional part that is made of the same material as an implant and abutment. It is applicable as an abutment for immediate loading, transfer for open/close tray, carrier for implant insertion, or healing abutment.

For multi-unit cases use M1+PCOM or M1 Multi-Unit and for a single crown choose from PCO titanium base.

Pre-milled abutment PMAB is a customizable one-piece abutment for a single crown metal or plastic PEEK framework.

ROOTT R abutment characteristics

- Angled access for tunnel from 0° to 20°;
- For bridges (PCOR, MS1, M1, CRE,) and single crowns (PCO1-PCO3, PMAB);
- Variety of gingiva part height PCO1-PCO3S;
- Variety of titanium bases height PCO and PCOS.

PMAB abutment characteristics

- Approved for use with a MEDENTIKA PreFace® Abutment Blank Holders;
- Provides unlimited possibilities to create high precision one-piece customized titanium abutment with an in house milling machine;
- Ideal adjustment for shape, emergence profile, esthetic properties - are available for frequently situation.




Diagram illustrating the Root R implant system, showing four variations (PCOR, PCO1, PCO2, PCO3) with increasing gingiva part heights (1.5mm, 2.5mm, 3.5mm). The diagram includes dimensions for the PCOR implant: 6mm height, 4mm diameter, 21mm length, and 29mm diameter. The legend indicates that the numbers mean gingiva part height.

REF	PCOR	PCO1	PCO2	PCO3
H	1.5 mm	1.5 mm	2.5 mm	3.5 mm

Instruction for Exocad using ROOTT R

Extraoral scan using PCO, PCOS

Step 1. Upload .slt file to Exocad software. Extraoral scan posts shall be visible.





Step 2. Choose abutment that will replace scan post.


Implant + abutment, material

e.g. ROOTT R + PCO1 zro

TRATE® ROOTT® - ROOTT R - PCO1 zro



**Intraoral or extraoral scan,
screw channel angle**
e.g. PCO1 25µm Extraoral
scan 0°



Cement gap size in microns, screw channel

e.g. PCO1 25µm screw
channel standard.

Step 3. Mark an area to detect Scan post position. Detected scan post shall change a color.

TRATE® ROOTT® - ROOTT R - MS1 PCOMS peek pmma

TRATE® ROOTT® - ROOTT R - PCOMS MS1 zro

TRATE® ROOTT® - ROOTT R - PCO1S peek pmma


TRATE® ROOTT® - ROOTT R - PCO1S zro

TRATE® ROOTT® - ROOTT R - PCO1 peek pmma

TRATE® ROOTT® - ROOTT R - PCO1 zro

TRATE® ROOTT® - ROOTT R - PCO2S peek pmma

TRATE® ROOTT® - ROOTT R - PCO2S zro



ROOTT R

Material – PEEK, PMMA or ZrO

Wizard Detect Implant Position ?

Tooth 36

SELECT IMPLANT

Full library


Favorites only

ROOTT R - PCO1 zro

PCO1 25µm Scr


PCO1 25µm Extr

Info about this library...



Click 'Best fit matching' when done


Best fit matching



Necessary products to make a prosthesis



Abutment
PCO



Screw
SFPCO




Screwdriver
SDLB




Extraoral scan using PMAB

Step 1. Upload .stl file to Exocad software. Extraoral scan post SPCO shall be visible (p. 14, Step 1.)

Step 2. Choose Pre-milled abutment PMAB that will replace scan post SPCO.



Step 3. Create a framework and mill the abutment.




Necessary products to make a prosthesis



Intraoral scan

Step 1. Upload scanned model to Exocad software. Intraoral scan post shall be visible.



Step 2. Choose abutment that will replace scan post (see p. 11, Step 2.).

Implant + abutment, material

e.g. ROOTT R + MS+ PCOMS zro

TRATE® ROOTT® - ROOTT R - MS1 PCOMS zro


SELECT IMPLANT

Full library Favorites only

TRATE® ROOTT® - ROOTT R - ▼ ★

PCOMS 35µm S▼ PCOMS 35µm Ir▼

Info about this library...



Intraoral or extraoral scan, screw channel angle

e.g. PCOMS 35 µm SPCOMIOS MS1 scan 0°

- PCOMS MS1 35µm SPCOMS MS1 scan 10°
- PCOMS MS1 35µm SPCOMS MS1 scan 15°
- PCOMS MS1 35µm SPCOMS MS1 scan 20°
- PCOMS MS1 35µm SPCOMIOS MS1 scan 0°**
- PCOMS MS1 35µm SPCOMIOS MS1 scan 5°
- PCOMS MS1 35µm SPCOMIOS MS1 scan 10°
- PCOMS MS1 35µm SPCOMIOS MS1 scan 15°
- PCOMS MS1 35µm SPCOMIOS MS1 scan 20°

PCOMS 25µm Screw channel standard

PCOMS 35µm Screw channel standard

PCOMS 50µm Screw channel standard

PCOMS 25µm Screw channel narrow

PCOMS 35µm Screw channel narrow


PCOMS 50µm Screw channel narrow

Cement gap size in microns, screw channel size

e.g. PCOMS 35 µm Screw channel standard

Intraoral scan

Step 1. Upload scanned model to Exocad software. Intraoral scan post shall be visible.



Step 2. Choose abutment that will replace scan post (see p. 11, Step 2.).

Implant + abutment, material

e.g. ROOTT R + MS+ PCOMS zro

TRATE® ROOTT® - ROOTT R - MS1 PCOMS zro


SELECT IMPLANT

Full library Favorites only

TRATE® ROOTT® - ROOTT R - ▼ ★

PCOMS 35µm S▼ PCOMS 35µm Ir▼

Info about this library...



Intraoral or extraoral scan, screw channel angle

e.g. PCOMS 35 µm SPCOMIOS MS1 scan 0°

- PCOMS MS1 35µm SPCOMS MS1 scan 10°
- PCOMS MS1 35µm SPCOMS MS1 scan 15°
- PCOMS MS1 35µm SPCOMS MS1 scan 20°
- PCOMS MS1 35µm SPCOMIOS MS1 scan 0°**
- PCOMS MS1 35µm SPCOMIOS MS1 scan 5°
- PCOMS MS1 35µm SPCOMIOS MS1 scan 10°
- PCOMS MS1 35µm SPCOMIOS MS1 scan 15°
- PCOMS MS1 35µm SPCOMIOS MS1 scan 20°

PCOMS 25µm Screw channel standard

PCOMS 35µm Screw channel standard

PCOMS 50µm Screw channel standard

PCOMS 25µm Screw channel narrow


PCOMS 35µm Screw channel narrow

PCOMS 50µm Screw channel narrow


Cement gap size in microns, screw channel size

e.g. PCOMS 35 µm Screw channel standard

New digital Multi Unit solution with ROOTT R




Ištrinsiu arba
pakeisiu konceptą




ROOTT **M** ROOTT **P**

Scan posts



ROOTT M Scan posts corresponds with multiunit platform connection. Scan posts indicate the exact position of the implant in the jaw. During the scanning process, the information about the position is transferred into digital format. According to Scan post position, height and direction, Scan-post is converted to abutment and sets the precise location of analog. Get yourself acquainted with all types of ROOTT M Scan posts, which are developed for the precise manufacturing of bridge prosthesis.



SPCOM scan post characteristic

- Extraoral scanning
- Long body allows comfortable usage and precise results of laboratory workflow
- Easily scannable
- Reusable



SPCOMIO scan post characteristic

- Intraoral scanning
- Short body allows comfortable usage and precise results of implantologists' workflow
- Easily scannable
- Reusable

Abutments

ROOTT M and ROOTT P digital library have two options for abutments. There are an abutment PCOM and a digital abutment that can be used directly without a physical one.

ROOTT M ROOTT P abutment characteristics

- Angled access for tunnel from 0° to 20°;
- Only for bridges.

TRATE® ROOTT® - ROOTT M ▾ ★

TRATE® ROOTT® - ROOTT M P - Multi - Unit

TRATE® ROOTT® - ROOTT M - PCOM peek pmma

TRATE® ROOTT® - ROOTT M - PCOM zro

TRATE® ROOTT® - ROOTT R - CRE peek pmma

TRATE® ROOTT® - ROOTT R - CRE zro

TRATE® ROOTT® - ROOTT R - M1 PCOM peek pmma

TRATE® ROOTT® - ROOTT R - M1 PCOM peek pmma

TRATE® ROOTT® - ROOTT R - MS1 PCOMS peek pmma

TRATE® ROOTT® - ROOTT R - PCOMS MS1 zro



ROOTT M ROOTT P
Only in Digital library

TRATE® ROOTT® - ROOTT M - PCOM zro

TRATE® ROOTT® - ROOTT R - CRE peek pmma

TRATE® ROOTT® - ROOTT R - CRE zro

TRATE® ROOTT® - ROOTT R - M1 Multi - Unit


TRATE® ROOTT® - ROOTT R - M1 PCOM peek pmma

TRATE® ROOTT® - ROOTT R - M1 PCOM zro

TRATE® ROOTT® - ROOTT R - MS1 Multi - Unit

TRATE® ROOTT® - ROOTT R - MS1 PCOMS peek pmma

TRATE® ROOTT® - ROOTT R - MS1 PCOMS zro





ROOTT M ROOTT P
Material – peek, pmma or zro

Instruction for Exocad using ROOTT M ROOTT P

Extraoral scan

Step 1. Scanned model upload to Exocad software. Extraoral scan posts shall be visible.



Step 2. Choose abutment that will replace scan post.

Implant + abutment, material

e.g. ROOTT M + PCOM zro

TRATE® ROOTT® - ROOTT M - PCOM zro

PCOM 25µm Screw channel standard

PCOM 35µm Screw channel standard

PCOM 50µm Screw channel standard

PCOM 25µm Screw channel narrow

PCOM 35µm Screw channel narrow

PCOM 50µm Screw channel narrow

SELECT IMPLANT

Full library Favorites only

TRATE® ROOTT® - ROOTT M

PCOM 25µm Scr PCOM 25µm Ext

Info about this li

Click 'Best fit match'

Best fit n

Intraoral or extraoral scan, screw channel angle
e.g. PCOM 25 µm Extraoral scan 10°

PCOM 25µm Extraoral scan 0°

PCOM 25µm Extraoral scan 5°

PCOM 25µm Extraoral scan 10°

PCOM 25µm Intraoral scan 0°

PCOM 25µm Intraoral scan 5°


PCOM 25µm Intraoral scan 10°

PCOM 25µm Intraoral scan 15°


PCOM 25µm Intraoral scan 20°

PCOM 25µm Extraoral scan 15°


Cement gap size in microns, screw channel size
e.g. PCOM 25 µm screw channel standard
(Suitable for screw driver and screw)




Standard screw channels



Cement gap 25µm




Prosthesis in Exocad




Digital model with analog ANMD

Necessary products to make a prosthesis




Abutment
PCOM

+




Screw
SFPCOM

+




Screwdriver
SDLB

→




Created by Pesterev Evgeniy




Digital MU
abutment

+




Screw
SFPCOM

+



Screwdriver
SDLB

→



Created by Sokratis Agapiou

ROOTT S

Scan post



ROOTT S Scan-posts corresponds with small multiunit platform connection. Scan-posts indicate the exact position of the implant in the jaw. During the scanning process, the information about the position is transferred into digital format. According to Scan-post position, height and direction, Scan-post is converted to abutment and sets the precise location of analog. Get yourself acquainted with all types of ROOTT S Scan-posts, which are developed for the precise manufacturing of bridge prosthesis.

*screw-retained restrations with a wide and secure fixation screw



SPCOMS scan post characteristic

- Extraoral scanning
- Long body allows comfortable usage and precise results of laboratory workflow
- Easily scannable
- Reusable



SPCOMIOS scan post characteristic



- Intraoral scanning
- Short body allows comfortable usage and precise results of implantologists' workflow
- Easily scannable
- Reusable

Abutments

ROOTT S digital library has two option for abutments. There are an abutment PCOMS and a digital abutment that can be used directly without a physical one.


ROOTT S abutments characteristics

- Angled access for tunnel from 0° to 20°;
- Only for bridges.



ROOTT S

Material – peek, pmma or zro




ROOTT S

Material – peek, pmma or zro

Instruction for Exocad using ROOTT S

Extraoral scan

Step 1. Scanned model upload to Exocad software. Extraoral scan posts shall be visible.



Step 2. Choose abutment that will replace scan post.

Implant + abutment, material

e.g. ROOTT S + PCOMS peek pmma

TRATE® ROOTT® - ROOTT S - PCOMS peek pmma

Cement gap size in microns, screw channel

e.g. PCOMS 90 µm Screw channel standard

PCOMS 70µm Screw channel standard




PCOMS 90µm Screw channel standard

PCOMS 110µm Screw channel standard


PCOMS 70µm Screw channel narrow

PCOMS 90µm Screw channel narrow

PCOMS110µm Screw channel narrow

TRATE® ROOTT® - ROOTT S - 
PCOMS 90µm S  PCOMS 90µm E 

[Info about this library...](#)



Click 'Best fit matching' when done

Best fit matching

Intraoral or extraoral scan, screw channel angle.

e.g. PCOMS 90 µm Extraoral scan 0°

PCOMS 90µm Extraoral scan 0°

PCOMS 90µm Extraoral scan 5°

PCOMS 90µm Extraoral scan 10°

PCOMS 90µm Intraoral scan 0°



PCOMS 90µm Intraoral scan 5°

PCOMS 90µm Intraoral scan 10°


PCOMS 90µm Intraoral scan 15°

PCOMS 90µm Intraoral scan 20°


PCOMS 90µm Extraoral scan 15°

Step 4. Mark an area to detect Scan post position.**Intraoral scan****Step 1.** Upload scanned model to Exocad software. Intraoral scan post shall be visible.**Step 2.** Choose abutment that will replace scan post (see p. 7 , Step 2.).

Necessary products to make a prosthesis



ROOTT M P S digital product scheme



* MU abutment is only accessible in digital library with angulation option and used with SFPCOMS screw for ROOTT S, SFPCOM screw for ROOTT M/P.

ROOTT C ROOT CS ROOT B ROOT BS

Scanable superstructures



ROOTT C, CS, B, BS Scan-post give a wide range of options. There are 45 different ways to scan with intraoral and extraoral scanners – choose from TRA, HE TOEA, TOE, TOES, telescopic abutments or External platform varieties. Scan-posts have a few height options that open possibilities for different clinical cases and patients mouth


Transfers can be used as scan-posts that make workflow more precise and effective. Scan-posts indicate the exact position and the depth of the screwed implant in the jaw. During the scanning process, the information about the position is transferred into digital format. According to Scan-post position, Scan-post is converted to telescopic abutment and sets the precise location of analog. Get yourself acquainted with one-piece abutments which are suitable even for complex clinical cases.




TRA




HE



TOEA



TOE



TOES

External platform
7, 5, 4, 3 mmTCE, PCE
TCES, PCES
TCEXS, PCEXS

Abutments

ROOTT C, CS, B, BS abutments are metal and plastic connectors between an implant and a crown. They do not require a screw making it the best solution for time-saving and comfort. Both TCE and PCE are the same size, however, the material is different – titanium or peek. Plastic peek gives more amortization, cushioning in the tooth, which makes the bite more comfortable and reduces the risk of fracture.



TCE0
0 mm

TCE1
1 mm

TCE2
2 mm

TCE3
3 mm



TCEX0
0 mm

TCEX1
1 mm

TCEXS1
1 mm

TCEXS2
2 mm

Instruction for Exocad using ROOTT C CS B BS

Extraoral scan

Step 1. Scanned model upload to Exocad software. Extraoral scan posts shall be visible.



Step 2. Choose abutment that will replace scan post.

Abutment, material

e.g. External platform TCE1 PCE1 peek pmma

External platform - TCE1 PCE1 peek pmma

Cement gap size in microns, screw channel
e.g. TCE PCE 70 µm screw channel standard
(Another possibility- no screw inside the crown)

TCE1 PCE1 70µm


TCE1 PCE1 90µm

TCE1 PCE1 110µm

TCE1 PCE1 70µm no screw channel

TCE1 PCE1 90µm no screw channel

TCE1 PCE1 110µm no screw channel



Scan post

e.g. TCE0 PCE0

TCE1 PCE1 90µm TOES Intraoral scan

TCE1 PCE1 90µm ANE Intraoral scan

TCE0 PCE0

TCE1 PCE1


TCE2 PCE2

TCE3 PCE3



TCEXS1 PCEXS1

TCEXS2 PCEXS2










Step 3. Mark an area to detect Scan post position. Detected scan post shall change a color.












Step 4. Create a crown's anatomy, print a model and insert a digital analog in a digital model.



ROOTT^M ROOTT^P

 <p>Digital abutment MU</p> <p>ROOTT-M-P_TRATE_ROOTT_Multi-Unit</p>	 <p>Analog ANMD</p>  <p>Scan posts SPCOM SPCOMIO</p>  <p>Screw SFPCOM</p>  <p>Connection geometry *</p> <p>PCOM.03_BRIDGE_0.025_15 PCOM.03_BRIDGE_0.025_20 ...</p>
 <p>Abutment PCOM</p> <p>ROOTT-M-P_TRATE_ROOTT_PCOM</p>	 <p>Analog ANMD</p>  <p>Scan posts SPCOM SPCOMIO</p>  <p>Screw SFPCOM</p> <p>Connection geometry *</p> <p>PCOM.03_0.070_0bp PCOM.03_0.070_0sp PCOM.03_0.070_5bp ... PCOM.03_0.070_20 PCOM.03_0.090_0 ...</p>

ROOTT^S


 <p>Digital abutment MU</p> <p>ROOTT-S_TRATE_ROOTT_Multi-Unit</p>	 <p>Analog ANMS</p>  <p>Scan posts SPCOMS SPCOMIOS</p>  <p>Screw SFPCOMS</p>  <p>Connection geometry *</p> <p>PCOMS.03_BRIDGE_0.015_15 PCOMS.03_BRIDGE_0.015_20 ...</p>
 <p>Abutment PCOMS</p> <p>ROOTT-S_TRATE_ROOTT_PCOMS</p>	 <p>Analog ANMS</p>  <p>Scan posts SPCOMS SPCOMIOS</p>  <p>Screw SFPCOMS</p> <p>Connection geometry *</p> <p>PCOMS.03_0.070_0bp PCOMS.03_0.070_0sp PCOMS.03_0.070_5bp ... PCOMS.03_0.070_20 PCOMS.03_0.090_0 ...</p>

	Model	Material	Superstructure	Connection geometry	Abutment
	TCE0 PCE0	Analog ANE	HE-END SPCOIOE TOE TOEA TOES TRA External platform 7, 5, 4, 3 mm TCE0 TCE1 TCE2 TCE3 TCES0 TCES1 TCES2 TCXS1 TCXS2	TCE0.01_0.070 TCE0.01_0.090	TCE0
	TCE1 PCE1	Analog ANE	HE-END SPCOIOE TOE TOEA TOES TRA External platform 7, 5, 4, 3 mm TCE0 TCE1 TCE2 TCXS1 TCXS2	TCE1.01_0.070 TCE1.01_0.090	TCE1
	TCE2 PCE2	Analog ANE	HE-END SPCOIOE TOE TOEA TOES TRA External platform 7, 5, 4, 3 mm TCE0 TCE1 TCE2 TCXS1 TCXS2	TCE2.01_0.070 TCE2.01_0.090	TCE2
	TCE3 PCE3	Analog ANE	HE-END SPCOIOE TOE TOEA TOES TRA External platform 7, 5, 4, 3 mm TCE0 TCE1 TCE2 TCXS1 TCXS2	TCE3.01_0.070 TCE3.01_0.090	TCE3
	TCES0 PCES0	Analog ANE	HE-END SPCOIOE TOE TOEA TOES TRA External platform 7, 5, 4, 3 mm TCE0 TCE1 TCE2 TCXS1 TCXS2	TCES0.01_0.070 TCES0.01_0.090	TCES0
	TCES1 PCES1	Analog ANE	HE-END SPCOIOE TOE TOEA TOES TRA External platform 7, 5, 4, 3 mm TCE0 TCE1 TCE2 TCXS1 TCXS2	TCES1.01_0.070 TCES1.01_0.090	TCES1
	TCES2 PCES2	Analog ANE	HE-END SPCOIOE TOE TOEA TOES TRA External platform 7, 5, 4, 3 mm TCE0 TCE1 TCE2 TCXS1 TCXS2	TCES2.01_0.070 TCES2.01_0.090	TCES2
	TCXS1 PCEXS1	Analog ANE	HE-END SPCOIOE TOE TOEA TOES TRA External platform 7, 5, 4, 3 mm TCE0 TCE1 TCE2 TCE3 TCES0 TCES1 TCES2 TCXS1 TCXS2	TCXS1.01_0.070 TCXS1.01_0.090	TCXS1
	TCXS2 PCEXS2	Analog ANE	HE-END SPCOIOE TOE TOEA TOES TRA External platform 7, 5, 4, 3 mm TCE0 TCE1 TCE2 TCE3 TCES0 TCES1 TCES2 TCXS1 TCXS2	TCXS2.01_0.070 TCXS2.01_0.090	TCXS2


Model Creator

An accurate printed model is needed to make a precise prosthesis. Therefore, it is crucial to know Model Creator software.

Note 1. Lift up a model in a way that analog fits in it. If analogs are visible, jaw model should be higher. If analogs are invisible, model is in the right position.



Incorrect- model is too low



Correct- good model position